



# Recombinant *Xenopus tropicalis* Transcription factor Sox-9 (sox9)

<b>Product Code</b>	CSB-BP738696XBF
<b>Abbreviation</b>	sox9
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6F2E7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Xenopus tropicalis</i> (Western clawed frog) ( <i>Silurana tropicalis</i> )
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNLLDPFMKM TEEQDKCMMSG APSPTMSEDS AGSPCPSGSG SDTENTRPQE NTFPKGDPPEL KKETEDEKFP VCIREAVSQV LKGYDWTLVP MPVRVNGSSK SKPHVKRPMN AFMVWAQAAR RKLADQYPHL HNAELSKTLG KLWROLLNEGE KRPFVEEAER LRIQHKKDHP DYKYQPRRRK SVKNGQSEQE DGAEQTHISP NAIFKALQAD SPSASSMSE VHSPGEHSGQ SQGPPTPPTT PKTDVQPGKP DLKREGRPLQ ESGRQPPHID FRDVDIGELS SEVISTIETF DVNEFDQYLP PNGHPGVGST QAPYTGSYGI NSTPSATPGA GPAWMSKQQQ QQQQQPQPPQ HSLSTINSEQ SQQSQRTHIK TEQLSPSHYS DQQQQHSPQQ LNYSSFNLQH YSSSYPTITR AQYDYTEHQG SNSYYSHASG QNSGLYSNFS YMNPSQRPMY TPIADTTGVP SIPQTHSPQH WEQPVYTQLT RP
<b>Source</b>	Baculovirus
<b>Target Names</b>	sox9
<b>Protein Names</b>	Recommended name: Transcription factor Sox-9 Short name= Xt-sox9
<b>Expression Region</b>	1-482
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.