



Recombinant *Candida glabrata* RNA polymerase II transcription factor B subunit 3 (TFB3)

Product Code	CSB-MP739452CZI
Abbreviation	TFB3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6FMP4
Product Type	Recombinant Protein
Immunogen Species	<i>Candida glabrata</i> (strain ATCC 2001 / CBS 138 / JCM 3761 / NBRC 0622 / NRRL Y-65) (Yeast) (<i>Torulopsis glabrata</i>)
Purity	>85% (SDS-PAGE)
Sequence	MSSAVASKGM NQGSEFDDNK DMCPICKTDR YLSPDVRFLV NPECYHKICE SCVDRIFFSLG PAPCPYKSCD KILRKNKFKT QIFDDVGVEK EVDIRKRVFN VFNKTLEDFD NDLEAYNKYL EEVEDIVYKL DNKIDVVETE EKLRTYEELN KQLIMNNMER SKKDLENFEQ RQQFEKEMRM KKRMLERQIE SEDKMNKEWA KREIVNRLAS STNADDVIEG VKNTVKLKKS SARRKLEEIN RVLQNNPYSS LMLSQGLSRK EDNVPFTPFN GDRDLDKRYT IDQESYDDPF IKDLENRKDY IASGFRSDFV YDRMLTEAFM GLGCVVAEEV
Source	Mammalian cell
Target Names	TFB3
Protein Names	Recommended name: RNA polymerase II transcription factor B subunit 3 Alternative name(s): RNA polymerase II transcription factor B 38 kDa subunit RNA polymerase II transcription factor B p38 subunit
Expression Region	1-330
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.