



Recombinant Bartonella henselae 50S ribosomal protein L15 (rpLO)

Product Code	CSB-MP739837BSG
Abbreviation	rpLO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6G2Y4
Product Type	Recombinant Protein
Immunogen Species	Bartonella henselae (strain ATCC 49882 / DSM 28221 / Houston 1) (Rochalimaea henselae)
Purity	≥85% (SDS-PAGE)
Sequence	MKLNELSDCK GAIRNRKRVG RGIGSGTGKT GGRGVKQQRSGVSLNGFEGGQMPIYRRL PKRGFRNFFA KDYNEVSLGR IQSAIDTGKL NIEKPVDIIA LKEAGIIRRV KDGVRLLSDG ELKVKISFNV SYASKAARLK IEKAGGQVHF PETV
Source	Mammalian cell
Target Names	rpLO
Protein Names	Recommended name: 50S ribosomal protein L15
Expression Region	1-154
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.