



# Recombinant *Xenopus laevis* Cyclin-L1 (ccn1)

<b>Product Code</b>	CSB-MP740099XBE
<b>Abbreviation</b>	ccn1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6GN15
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Xenopus laevis</i> (African clawed frog)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAAVPQLSAP SAPARSADGI LIGDRQYSEV YLTIDYSLIP EERLSPTPSM SDGLDLNTET DLRILGCELI QSAGILLRLP QVAMATGQVL FHRFFYSKSF VKHSFEIIAM ACINLASKIE EAPRRIRDVI NVCHHLRQIR AKRTPSPLIL DQSYINTKNH VIKAERRILK ELGFCVHVKH PHKIIVMYLQ VLECERNQTL VQTAWNYMND CLRTNVFVRF DAETIACACI YLAARALQLS LPNRPHWFL FGATEENIQD ICITTLRLYS RIKPNYEFLE KEVDKRKVAL QEAKLKAKGL NPDGTPAILS MGGFSPASKP SSPRDVKTEE KSPNFAKVKR EMDDKQSSKS PYNGLRKENK RSRSVSRSR SRTKSRSRSHS PRRHYNRRR SRSGTYSSRS RSRSRSHSES PRRHHNHGSP HMKLKHRVED LRGRHAHKRK KSHSPSKSRE PSELAKKHRH EHGHRERRE RSRSFERSHK NKHHGSSHSG HGRHRR
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ccn1
<b>Protein Names</b>	Recommended name: Cyclin-L1
<b>Expression Region</b>	1-496
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.