



Recombinant Rat Keratin, type I cytoskeletal 39 (Krt39)

Product Code	CSB-BP740232RA
Abbreviation	Krt39
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6IFW3
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	MDTKGSTVTI SSSTPPQNCS GNTNVRTNSS NKSCYHDVQS TGHALQTPQG QGCRPSPCLY RCPNYLIRTY SFHPCPDDCS RCSDGINSHE KETMQILNER LASYLEKVRM LEGENADLED KIQEECSKTL PILCPDYSY YTTIEQLQK ILCTKAENSR LVSQIDNTKL AADDLRAKYE AELSLRQLVE ADANGLKQIL DALTLKADL EARVQSLTEE LLCLKTNHEE EINSLQCQLG DRINIEVTAA PSVDLNQILQ KMRCQYESIV ETNRKDVEEW FNTQMEELNQ QVVSSSQQQ CCQKDIIELR RTISALEVEL QAQHRMRDSQ ECILAETEAR YTALLAQIQS LIHNLEAQVA EIRSALQRQN QEYEVLLDIK SRLECEIATY RSLLESLDGR LPCNPCTTTW EPSCQARAME CLTPVYTSIS LPGIHKPCRA SGPPSRILVK ICTITKEIKD GKVISSHEHV QPCYITRPAK V
Source	Baculovirus
Target Names	Krt39
Protein Names	Recommended name: Keratin, type I cytoskeletal 39 Alternative name(s): Cytokeratin-39 Short name= CK-39 Keratin-39 Short name= K39 Type I hair keratin Ka35
Expression Region	1-481
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a member of the type I (acidic) keratin family, which belongs to the superfamily of intermediate filament (IF) proteins. Keratins are heteropolymeric structural proteins which form the intermediate filament. These filaments, along with actin microfilaments and microtubules, compose the cytoskeleton of epithelial cells. The type I keratin genes are clustered in a region of chromosome 17q12-q21.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.