



Recombinant Mouse Keratin, type I cytoskeletal 40 (Krt40)

Product Code	CSB-MP740233MO
Abbreviation	Krt40
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6IFX3
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MASDGSPSCC SSEPCAGASG CATASFCPTN TTCLPNTCST SRCQTSPFLC RASLPACCLS PCYLAGGCNS PCLVGSCAWC EEGSFNSNEK ETMQFLNDRL ASYLERVRSLEENNAELECR IREQCEPNAP LVCPDYQRYF DTIEELQQKI LCTKAENSRL AVQVDNCKLA ADDFRSKYES ELSLRQLVEN DITGLRGILG ELTLCKSDLE AHVESMKDDL ICLKKGHEEE VNULLREQLGD RLSVELDTAP TVDLNKVLDE MRCQYERVLA NNRRDAEEWF AAQTEELNQQ QKSSAEQLEG CQTEMLELKR KANTLEIELQ AQQTLTESLE CTVAETEAQY STQLAQMQL IDSVEHQLAE IRCDLERQNG EYQVLLDTKA RLECEINTYR GLLEKEDSRL PCNPGSGAPM PNSTCEPCSN SMCEPCSAYV ICTVENCCA
Source	Mammalian cell
Target Names	Krt40
Protein Names	Recommended name: Keratin, type I cytoskeletal 40 Alternative name(s): Cytokeratin-40 Short name= CK-40 Keratin-40 Short name= K40 Type I hair keratin Ka36
Expression Region	1-439
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a member of the type I (acidic) keratin family, which belongs to the superfamily of intermediate filament (IF) proteins. Keratins are heteropolymeric structural proteins which form the intermediate filament. These filaments, along with actin microfilaments and microtubules, compose the cytoskeleton of epithelial cells. The type I keratin genes are clustered in a region of chromosome 17q12-q21.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.