



Recombinant Photobacterium profundum Lipoyl synthase (lipA)

Product Code	CSB-MP740382PIG
Abbreviation	lipA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6LN87
Product Type	Recombinant Protein
Immunogen Species	Photobacterium profundum (strain SS9)
Purity	>85% (SDS-PAGE)
Sequence	MSKPIKIEQG IKYRDADKMA LIPVRNVAEE APKEVLRKPA WMKIKLPSSD KRIQEIKSAL RKNKLHSVCE EASCPNLAEC FNHGTATFMI LGAICTRRCP FCDVAHGRPL PPNAEESH L AQTADMKLK YVITSVDRD DLRDGGGAQHF VDCIREIREK SPEIHITLV PDFRGRMDRA LDILQGTPPN VFNHNLETAP RLYRKARPGA NYQWSLDLLK NFKEIHPEVP TKSGVMMGLG ETKEEIIQVL KDLRAHGVTM LTLGQYLAPS RHHLPVERYV PPAEFDELKE IALELGFTHA ACGPFVRSY HADLQAQGLE IK
Source	Mammalian cell
Target Names	lipA
Protein Names	Recommended name: Lipoyl synthase EC= 2.8.1.8 Alternative name(s): Lip-syn Short name= LS Lipoate synthase Lipoic acid synthase Sulfur insertion protein LipA
Expression Region	1-322
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.