



Recombinant Rat Tripartite motif-containing protein 40 (Trim40)

Product Code	CSB-MP740497RA
Abbreviation	Trim40
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6MFY8
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MVPLDKDNQD ICPICLDPLK EAVSTDCRHL FCRMCLIRHM DKASVSGVLS CPVCRKPCSE TVLGDNYICH THQKRVCRFC ESSRHLLCEE CLQSPEHRAH TELSIENAIH HYKERLNRRS RKLKRD LGDL QRLKAQEEKM LQALQVDWGS HRPRTEQQNQ DQTELQQKAL PRHWLDQRED PPEEVAKVFN FSEAVTQLSI LVSSLERMAK ELDASTLKDA SDLLDRSSQQ KLEGLLSHVP PANSKLS
Source	Mammalian cell
Target Names	Trim40
Protein Names	Recommended name: Tripartite motif-containing protein 40 Alternative name(s): RING finger protein 35
Expression Region	1-247
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.