



Recombinant *Bdellovibrio bacteriovorus* DNA-directed RNA polymerase subunit alpha (rpoA)

Product Code	CSB-BP740513BFP
Abbreviation	rpoA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6MJ36
Product Type	Recombinant Protein
Immunogen Species	<i>Bdellovibrio bacteriovorus</i> (strain ATCC 15356 / DSM 50701 / NCIB 9529 / HD100)
Purity	>85% (SDS-PAGE)
Sequence	MQEHYKFWR EMIKPKGFEV DRDSL RDDYA KFIIRPLERG FGVTLGNSLR RILLSSMMGS AITAVKFEGV LHEFTTIPDV LEDVTDIILN LKEVRFKQYT ADSLTLKISK KGPGKVTAAD IQTTDKIEVL NPDLP IATLG ANANFNAEIV VSFGRGYVPV ENRETDLPVG FIGVDALYSP IRKVNYNVS N ARVGQR TDYD ALTLEVWTDG SLKPEEAVAL SSKIMKEQLQ IFLTFDETME PAEEARELGS PTLNENLFRS VDDLELSVRS ANCLKNANIR YIGELVVRSE AEMLKTKNFG RKSLNEIKEI LGEMGLGLGM KIEGWPPAGW DPSQPPQKRE TQQ
Source	Baculovirus
Target Names	rpoA
Protein Names	Recommended name: DNA-directed RNA polymerase subunit alpha Short name= RNAP subunit alpha EC= 2.7.7.6 Alternative name(s): RNA polymerase subunit alpha Transcriptase subunit alpha
Expression Region	1-343
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.