



Recombinant Arabidopsis thaliana Tubby-like F-box protein 11 (TULP11)

Product Code	CSB-MP740640DOA
Abbreviation	TULP11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NPQ1
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	<p> MTLRSLILEM RSRPHRVVHD LAAAAAADST SVSSQDYRWS EIPEELLREI LIRVEAADGG GWPSRRSVVA CAGVCRGWRL LMNETVVVPE ISSKLTFPIS LKQPGPRDSL VQCFIKRNRI TQSYHLYLGL TNSLTDDGKF LLAACKLKHT TCTDYIISLR SDDMSRRSQA YVGKVRSNFL GTKFTVFDGN LLPSTGAAKL RKRSYNPAK VSAKVPLGSY PVAHITYELN VLGSRGPRKM QCLMDTIPTS TMEPQGVASE PSEFPLLGTR STLSRSQSKP LRSSSSHLKE TPLVLSNKTP RWHEQLRCWC LNFHGRVTVA SVKNFQLVAA GASCGSGTGM SPERQSERII LQFGKVGKDM FTMDYGYPIA AFQAFACLS SFETRIACE </p>
Source	Mammalian cell
Target Names	TULP11
Protein Names	Recommended name: Tubby-like F-box protein 11 Short name= AtTLP11
Expression Region	1-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.