



Recombinant *Xenopus laevis* Polyribonucleotide 5'-hydroxyl-kinase Clp1 (clp1)

Product Code	CSB-BP740650XBE
Abbreviation	clp1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NS21
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	MSEDLSAPG PSAPDSSWSS TGPSTKFELE RETELRLEVE GTDPVRVELV SGLAEVFGTE LTRNKKYTFP PGSRAAIFTW HGCTVQLWGS PDMAYVSRDT PMLLYLNTQV GLEQMRVQAE REGERGPRVL VAGPSDVGKS TLCRLLNLYA VRRGRRPTLV ELDVGQGSVS VPGTMGALCV ERPADVEEGF SAQAPLVYHF GSTTPGTNIK LYNKLTSRLA HVFNLRCDSN RRASVSGCLI NTCGWVKGSG YQALIIHAASA FEVDVVLVLD QERLYNDLLR DLPHFVRTLL LPKSGGASER SKECRRESRD QRVREYFYGP RGSLYPHAFE IKFSEVRVYK VGAPTIPDSC LPLGMSQEDN QLKLVPVTPG RDMAHHLLSV VPLDGGSAEE GIEERSVAGF IVITGV DTER QTLTLLSPAP RPLPKCVLLI MDIRFMDLK
Source	Baculovirus
Target Names	clp1
Protein Names	Recommended name: Polyribonucleotide 5'-hydroxyl-kinase Clp1 EC= 2.7.1.78 Alternative name(s): Polynucleotide kinase Clp1 Pre-mRNA cleavage complex II protein Clp1
Expression Region	1-439
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.