



# Recombinant Human EF-hand domain-containing protein C3orf25 (C3orf25)

<b>Product Code</b>	CSB-BP740681HU
<b>Abbreviation</b>	C3orf25
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6NXP0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDDDYEAYHS LFLSLLGLCP SKTPINENAP VFDPEPVIAH CFKQFQQKDF RLPQTRRRRII MVPRKEDQTP LNPASQPQAP PKPIPSFKVL EARDIQEQPE DRKTWLSQRS KLRQELESFG DVKRWLENKP SITPSEAKVL HMIHEEQSAQ PNASQATTRT TRKKAPRLSR LSRQMVPQLQ LPEPPALSVM YSYLHSRKIK ILEIFHKVGQ GENQRITREE FIAAVKAVGV PLKNQEVEDI VIYLSSLGKH NTITMDILAN TYKQWSMAQQ RSSLATAREH YILAKHRDSL KGPLKKQEVD SAPQLPKVDL LTVPAVDTQM ETRPMTLEEM EEVGKRYRER QRQHKLTIPS IQYTEQCHLV RCGNRHFDEH CLPSTIHGDM RELIDSARRH NFLVYLQCWK LCKSYGLPLT EDILMKALLY PGDKIIFQMD KVCPIRQPGG YYSDWKVFSP NLALLRSQGP GKSKRTDKKT PKKSKKMRFK EFEEFTRKLLK VKRSSGLQQT HPNSFWPGHL LDKLQLYLPT VATDRSLALF SCVQHQP HVY PATYHPDHHWW PLRNKNYMTH AHYDAAKVYY IN
<b>Source</b>	Baculovirus
<b>Target Names</b>	EFCAB12
<b>Protein Names</b>	Recommended name: EF-hand domain-containing protein C3orf25
<b>Expression Region</b>	1-572
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.