



# Recombinant Human Methyltransferase-like protein 2B (METTL2B)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-YP740703HU  |
| <b>Abbreviation</b>      | METTL2B   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q6P1Q9  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Homo sapiens (Human)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | AGSYPEGAP AILADKRQQF GSRFLSDPAR VFHHNAWDNV EWSEEQAAAA<br>ERKVQENSIQ RVCQEKQVDY EINAHKYWND FYKIHENGFF KDRHWLWLFTEF<br>PELAPSQNQN HLKDWFLENK SEVCECRNNE DGPGLIMEEQ HKCSSKSLEH<br>KTQTPPVEEN VTQKISDLEI CADEFPGSSA TYRILEVGCG VGNTVFPILQ<br>TNNDPGLFVY CCDFSSTAIE LVQTNSEYDP SRCFAFVHDL CDEEKSYVPV<br>KGSLDIIILI FVLSAVVPDK MQKAINRLSR LLKPGGMVLL RDYGRYDMAQ<br>LRFKKGQCLS GNFYVRGDGT RYVFFTQEEL DTLFTTAGLE KVQNLVDRRL<br>QVNRGKQLTM YRVWIQCKYC KPLLSSTS |
| <b>Source</b>            | Yeast   |
| <b>Target Names</b>      | METTL2B   |
| <b>Protein Names</b>     | Recommended name: Methyltransferase-like protein 2B EC= 2.1.1.-   |
| <b>Expression Region</b> | 2-378   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |
| <b>Target Details</b>    | This gene is a member of a family of methyltransferases that share homology with, but are distinct from, the UbiE family of methyltransferases. Alternatively spliced variants which encode different protein isoforms have been described; however, not all variants have been fully characterized.  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients,   |



storage temperature and the stability of the protein itself.  
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.