



Recombinant Human Serine/threonine-protein kinase Nek5 (NEK5)

Product Code	CSB-BP740713HU
Abbreviation	NEK5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6P3R8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDKYDVIKAI GQGAFGKAYL AKGKSDSKHC VIKEIFEKM PIQEKEASKK EVILLEKMKH PNIVAFFNSF QENGRLFIVM EYCDGGDLMK RINRQRGVLF SEDQILGWV FV QISLGLKHIH DRKILHRDIK AQNIFLSKNG MVAKLGDFGI ARVLNNSMEL ARTCIGTPYY LSPEICQNKP YNNKTDIWSL GCVLYELCTL KHPFEGNNLQ QLVLKICQAH FAPISPGFSR ELHSLISQLF QVSPRDRPSI NSILKRPFLE NLIPKYLTPE VIQEEFSHML ICRA GAPASR HAGKVVQKCK IQKVRFQGKC PPRSRSIVPI KR NAILHRNE WRPPAGAQA RSIKMIERP IAAVCGHYDY YYAQLDMLRR RAHKPSYHPI PQENTGVEDY GQETRHGPSP SQWPAEYLQR KFEAQYKYLK VEKQLGLRPS SAEPNYNQRQ ELRSNGEPR FQELPFRKNE MKEQEYWKQL EEIRQQYHND MKEIRKKMGR EPEENSKISH KTYLVKKS NL PVHQDASEGE APVQMEFRSC CPGWSAMARS WLTATSASQD IEKDLKQMLR QNTKESKNPE QKYKAKKGVK FEINLDCIS DENILQEEEA MDIPNETLTF EDGMKFKEYE CVKEHGDYTD KAFEKLHCPE AGFSTQTVAA VGNRRQWDGG APQTLQMMMA VADITSTCPT GPDSSEVLSV SRQEGKTKDP YSPVLILM
Source	Baculovirus
Target Names	NEK5
Protein Names	Recommended name: Serine/threonine-protein kinase Nek5 EC= 2.7.11.1 Alternative name(s): Never in mitosis A-related kinase 5 Short name= NimA-related protein kinase 5
Expression Region	1-708
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.