



Recombinant Danio rerio Pre-mRNA-splicing factor 38B (prpf38b), partial

Product Code	CSB-BP740740DIL
Abbreviation	prpf38b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6P7Y3
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MAGSQQQQQQ QAVSKPTGGK HGNNLPLWGN EKTMNLNPMI LTNVLSSPYF KVQLYELKTY HEVVDEIYFK VNHVEPWKEG SRKTAGQTGM CGGVRGVGTG GIVSTAFCLL YKLFCLKLTR KQVMGLITHS DSPDIRALGF MYIRYTPPPP DLVDWYDEFL DDEEELDVKA GGGCVMTVGE MLRSFLTKLE WFSTLFPRIIP VPVQKAIDQH MKSRPRKPPQ KDEQEEEEEE ATAAPAADTG RHGDKRRSRT PRRSPSPRKS QNRSRSRSHH RERHGASFDY ELERERDRQR KEREGKDRDR DRDRDRERDR ERDRDRRRSR TPDRNAERRR SRSRERRRRS STSRDKRTER KDRDKDREAE SERERSRKKD REHHKDRERS KDKRSKGEGE ERRHKDDKEE KKHREEKRSK RSRSRSRDRK HKAERSSKKR SRSGRSRQEQ AGEEKNRKRE RSHSKDRQHK RSRSKERSHR RESSNERIHA RQERPSSSESG ERTNSVRADS P
Source	Baculovirus
Target Names	prpf38b
Protein Names	Recommended name: Pre-mRNA-splicing factor 38B
Expression Region	1-501
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.