



# Recombinant Danio rerio Pre-mRNA-splicing factor 38B (prpf38b), partial

<b>Product Code</b>	CSB-MP740740DIL
<b>Abbreviation</b>	prpf38b
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6P7Y3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAGSQQQQQQ QAVSKPTGGK HGNNLPLWGN EKTMNLNPMI LTNVLSSPYF KVQLYELKTY HEVVDEIYFK VNHVEPWKEG SRKTAGQTGM CGGVRGVGTG GIVSTAFCLL YKLFCLKLTR KQVMGLITHS DSPDIRALGF MYIRYTPPPP DLVDWYDEFL DDEEELDVKA GGGCVMTVGE MLRSFLTKLE WFSTLFPRIIP VPVQKAIDQH MKSRPRKPPQ KDEQEEEEEE ATAAPAADTG RHGDKRRSRT PRRSPSPRKS QNRSRSRSHH RERHGASFDY ELERERDRQR KEREGKDRDR DRDRDRERDR ERDRDRRRSR TPDRNAERRR SRSRERRRSR STSRDKRTER KDRDKDREAE SERERSRKKD REHHKDRERS KDKRSKGEGE ERRHKDDKEE KKHREEKRSK RSRRSRDRK HKAERSSKKR SRSGRSRQEQ AGEEKNRKRE RSHSKDRQHK RSRSKERSHR RESSNERIHA RQERPSSSESG ERTNSVRADS P
<b>Source</b>	Mammalian cell
<b>Target Names</b>	prpf38b
<b>Protein Names</b>	Recommended name: Pre-mRNA-splicing factor 38B
<b>Expression Region</b>	1-501
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.