



Recombinant Mouse WAS/WASL-interacting protein family member 2 (Wipf2)

Product Code	CSB-MP740777MO
Abbreviation	Wipf2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6PEV3
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	<p> MPIPPPPPPP PGPPPPTFN QANTEQPKLS RDEQRNRGAL LQDICKGTKL KKVTNVNDRS APVIEKPRGS SGGYGPGAAA LQPKGGFLFQG GVPKLRPVGA KDASEAPAGK PALQVPSSRA AAPRPPGSAA SGRPHDDTDS NRASLPELPR MQRPSLPDLS RPNTASGTGM KHSSSAPPPP PPGRRANAPP TPLPLHSNKA QAYNREKPLP PTPGQRLHPG REGHPAPPPV KPPSPVNIR TGPSGQSLAP PPPPYRQPPG VPNGPSSPTN ESAPELPQRH NSLHRKTPGP VRGLAPPPPT SATPSLLSNR PPPPAREPPS RGAAPPPPPP MIRNGARDAP PPPPPYRMHG SEPPSRGKPP PPPSRTPAGP PPPPPPPLRN GHRDSITTVR SFLDDFESKY SFHPVEDFPA PEEYKHLQRV YPSKTNRAAR GAPPLPPILR </p>
Source	Mammalian cell
Target Names	Wipf2
Protein Names	Recommended name: WAS/WASL-interacting protein family member 2 Alternative name(s): WASP-interacting protein-related protein WIP-related protein
Expression Region	1-440
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a WASP interacting protein (WIP)-related protein. It has been shown that this protein has a role in the WASP-mediated organization of the actin cytoskeleton and that this protein is a potential link between the activated platelet-derived growth factor receptor and the actin polymerization machinery.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.