



Recombinant Arabidopsis thaliana AP2-like ethylene-responsive transcription factor SNZ (SNZ)

Product Code	CSB-MP740803DOA
Abbreviation	SNZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6PV67
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MLDLNLGILS THNEDEDCKV PTSIFIQEED SINPSNDNLS LITFGILKRN VEILPPPPPP PPPPPPSENE LSGPGNEWLD LSSMQRNKQE TLVMKKKSRR GPRSRSSHYP GVTFYRRTGR WESHIWDCGK QVYLGGFDTA YTAARAYDRA AIRFRGLQAD INFIVDDYKQ DIEKMKNLSK EEFVQSLRRA SASLARGGSK YKNTHMRNDH IHLFQNRGLN AAAAKCNEIR KMEGDIKLGA HSGKNEHNDL ELSLGISSSS KVRILEPADY YMGLNRSVTS LHGKPLPGYL PITEIKPLKT VVASSGFPI TMINPSSLSL SCFDP
Source	Mammalian cell
Target Names	SNZ
Protein Names	Recommended name: AP2-like ethylene-responsive transcription factor SNZ Alternative name(s): Protein SCHNARCHZAPFEN
Expression Region	1-325
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.