



Recombinant Dog Cyclin-dependent kinase inhibitor 1B (CDKN1B)

Product Code	CSB-MP740867DO
Abbreviation	CDKN1B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6SLL5
Product Type	Recombinant Protein
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Purity	>85% (SDS-PAGE)
Sequence	MSNVRVSNNGS PTSERRDAKQ AEYPKPSACR NLFPGVNHEE LTRDLEKHCI DMEEASQNKW NFDFQNHKPL EGKYEWQEVE KGSLPEFYR PPRPPKGACK VPAQESQDVS GTRQAGPLLQ SQANSEDTLH VDQKTDAPDS QTGLAEQCTG IRKRPATDDS SPQNKRANRT EENVSDGSPN AGSVEQTPKK PGLRRRQT
Source	Mammalian cell
Target Names	CDKN1B
Protein Names	Recommended name: Cyclin-dependent kinase inhibitor 1B Alternative name(s): Cyclin-dependent kinase inhibitor p27 p27Kip1
Expression Region	1-198
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a cyclin-dependent kinase inhibitor, which shares a limited similarity with CDK inhibitor CDKN1A/p21. The encoded protein binds to and prevents the activation of cyclin E-CDK2 or cyclin D-CDK4 complexes, and thus controls the cell cycle progression at G1. The degradation of this protein, which is triggered by its CDK dependent phosphorylation and subsequent ubiquitination by SCF complexes, is required for the cellular transition from quiescence to the proliferative state.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.