



Recombinant Mouse Kelch-like protein 6 (Klh6)

Product Code	CSB-EP740929MO
Abbreviation	Klh6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6V595
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MAGQRGACTM SDVVERSLEG PLALTTDKSS QKRGDLVEIL TGEKVKFDDT RLSLILQNGL ETLREENALT DVILCVDVQE FSCHRVLAA ASNYFRAMFC NDLKEKEYEKR IIIKGVDAET MHTLLVYTYT SKALITKQNV QRVLEAANLF QFLQLVDACA SFLTEALNPE NCIGILRLAD THSLDSLKMQ VQSYIIQNFV QILNSEEFLE LPVDTLYHIL RSDDL CVTEE AQVFETVMSW VRHKQSERLC LLPCVLENVR LPLLDPWYFV EMVEADPLIR QCPEVFPLLQ EARMYHLSGN EIISSERTKPR MHEFQSEVFM IIGGCTKDER FVAEVTCLDP LRRSRLEVAK LPLTEHELES ENKKWVEFAC VTLKNEVYIS GGETQHDVW KYNSSINKWI QIEYLNIGRW RHKMVVLGGK VYVLGGFDGL QRINNVETYD PFHNCWSEAA PLHVHVSSFA ATSHKKKLYV IGGGPNGLA TDKTQCYDPS TNKWVLKSAM PVEAKCINAV SFQDHIYVVG GAMRALYAYS PLEDSWCLVT QLSHERASCG IAPCANNKLYI TGGRDEKNEV IATVLCWDPE AQLTTEECVL PRGVSHHGSV TIRKSYTHIR RVVPGAVSV
Source	E.coli
Target Names	Klh6
Protein Names	Recommended name: Kelch-like protein 6
Expression Region	1-619
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.