



Recombinant Human E3 ubiquitin-protein ligase NHLRC1 (NHLRC1)

Product Code	CSB-MP740936HU
Abbreviation	NHLRC1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6VVB1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAAEASESGP ALHELMREAE ISLLECKVCF EKFGHRQRR PRNLSCGHVV CLACVAALAH PRTLALCPFP CRRACRGCDT SDCLPVLHLI ELLGSALRQS PAAHRAAPSA PGALTCHHTF GGWGTLVNPT GLALCPKTGR VVVVHDGRRR VKIFDSGGGC AHQFGEKGD A QDIRYPVDV TITNDCHVVV TDAGDRSIKV FDFGQIKLV IGGQFSLPWG VETTPQNGIV VTDAEAGSLH LLDVDFAEV LRRTERLQAH LCNPRGVAVS WLTGAIIVLE HPLALGTGVC STRVKVFSS MQLVGQVDTF GLSLYFPSKI TASAVTFDHQ GNVIVADTSG PAILCLGKPE EFPVPMKMT HGLSHPVALT FTKENSLVL DTASHSIKVY KVDWG
Source	Mammalian cell
Target Names	NHLRC1
Protein Names	Recommended name: E3 ubiquitin-protein ligase NHLRC1 EC= 6.3.2.- Alternative name(s): Malin NHL repeat-containing protein 1
Expression Region	1-395
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.