



Recombinant Human Cell cycle checkpoint control protein RAD9B (RAD9B)

Product Code	CSB-MP740946HU
Abbreviation	RAD9B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6WBX8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MLKCVMSGSQ VKVFGKAVQA LSRISDEFWL DPSKKGLALR CVNSSRSAYG CVLFSPVFFQ HYQWSALVKM SENELDTTLH LKCKLGMKSI LPIFRCLNSL ERNIEKCRIF TRSDKCKVVI QFFYRHGIKR THNICFQESQ PLQVIFDKNV CTNTLMIQPR LLADAILVFT SSQEEVTLAV TPLNFCLKSS NEESMDLSNA VHSEMFVGSDF EFDFFQIGMD TEITFCFKEL KGILTFSEAT HAPISYFDF PGKPLALSID DMLVEANFIL ATLADEQSRA SSPQSLCLSQ KRKRSDLIEK KAGKNVTGQA LECISKKAAP RRLYPKETLT NISALENCGS PAMKRVDGDV SEVSESSVSN TEEVPGSLCL RKFSCMFFGA VSSDQQEHFN HPFDSLARAS DSEEDMNNVC CRKEFNGSDA KYFCII
Source	Mammalian cell
Target Names	RAD9B
Protein Names	Recommended name: Cell cycle checkpoint control protein RAD9B Alternative name(s): DNA repair exonuclease rad9 homolog B Short name= hRAD9B
Expression Region	1-426
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.