



Recombinant *Oryza sativa* subsp. japonica PHD finger protein ALFIN-LIKE 9 (Os07g0608400, LOC_Os07g41740)

Product Code	CSB-MP741012OFG
Abbreviation	Os07g0608400, LOC_Os07g41740
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6YTY3
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MDAQYNPRTV EEVFRDFKGR RAGLVRALTA DVEDFFRQCD PEKENLCLYG FPNEHWEVNL PAEEVPELP EPALGINFAR DGMQEKDWLS MVAVHSDAWL LSVAFYFGAR FGFDKNDKR LFGMINDLPT IFEVVSQKSK AKPPSANNHS NSKSKSSNKT KSSEPRAKQP KPQPQPPVKN EGREEGGPD DEEGGGGGGG GGREEEHGET LCGACGESYG ADEFWICCDI CEKWFHGKCV KITPAKAEHI KQYKCPSCSG GNGGGGGVSG NGKRARPS
Source	Mammalian cell
Target Names	Os07g0608400
Protein Names	Recommended name: PHD finger protein ALFIN-LIKE 9
Expression Region	1-278
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.