



Recombinant *Oryza sativa* subsp. *japonica* Putative heat stress transcription factor B-4a (HSFB4A)

Product Code	CSB-YP741036OFG
Abbreviation	HSFB4A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6Z9R8
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MEWEESEAA RQKAAAASAS VVPAPFLTKT YQLVDDPATD HVVSWEDDDG GESASSFVWV RPPEFARDIL PNYFKHSNFS SFVRQLNTYG FRKVVPERWE FANEFFRKGE KQLLCEIHRK KSAATWPPF PPPPPFFAP RHFAAGAFFR HGDGMLHGRL GALVTTTERR HWFESAALPV APSSRLLSQL GPVIAPARRA AATPEEEALM QENHRLLRGN AALVQELAHM RKLSDIIF VQNHVRPVAP SPAAAAALHG LGVLRPPPAG GKAPASEVRG ASGRSATSSS SLTVAEDQPT LLALRLPRTT EKIINEVSGG NGGGSTKLFV VHLSSADEQT SSGASRKRSP PQEQPPTSPA PKRTLVEHS ELRLSIVSPP
Source	Yeast
Target Names	HSFB4A
Protein Names	Recommended name: Putative heat stress transcription factor B-4a Alternative name(s): Heat stress transcription factor 20 Short name= OsHsf-20
Expression Region	1-380
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.