



Recombinant Human Ras and Rab interactor-like protein (RINL)

Product Code	CSB-MP741063HU
Abbreviation	RINL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ZS11
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAQPEDKAPE VPTEGVRLVP PQVKNADRTP LGVLSTLEPL TRLQRTWGVW HVP ELDTQDA EALVGLWPLG SFLVTGRDPS QALVLRSGPL PGEVNTYQIQ KIPRGVSLES SNLCMPDLPH LLAFLSASRD VLPRTLPP PTLGPRDEHT DPVQIGRVQQ DTPGKVLSIV NQLYLETHRG WGREQTPQET EPEAAQRHDP APRNPAPHGV SWVKGPLSPE VDHPGPALAS LLEEEEDLE GKEEGREDDP EEEGPEDVLT IHVQSLVRAR SSVARQYRS LRVRIASDSG GPHGSGDPAT ELLQDVRHLL TDLQDHLAKD SYIRAVFGSR GPGLPKKDED PGPALETAVC QAVLAPLKPA LWTRLRTLRA PELRRLRRRQ TALRAGAGPP GAQGPGPEGQ SPAPALRSRI HERLAHLHAA CAPRRKVALL LEVCRDVYAG LARGENQDPL GADAFLPALT EELIWSPDIG DTQLDVEFLM ELLDPDELRG EAGYYLTTFW GALHHIAHYQ PETDRAPRGL SSEARASLHQ WHRRRTLHRK DHPRAQANLP FKEPWAEETV TGTSDN
Source	Mammalian cell
Target Names	RINL
Protein Names	Recommended name: Ras and Rab interactor-like protein
Expression Region	1-566
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.