



Recombinant Human LanC-like protein 3 (LANCL3)

Product Code	CSB-MP741072HU
Abbreviation	LANCL3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ZV70
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDTKRCFANR FDDYQGSLLA GQCEEAVAPL VTATIERILQ ELPPLGGGAE ARGATAGASA CQGGLYGGVA GVAYMLYHVS QSPLFATARE RYLRSKRLLI DACARAEWEG EPDADTRAAF LLGGAGVYAV ATLVYHALGR SDYVQPLGKF RALCAVCAPV SFLECGSDEL FVGRAGYLCA ALVLKQKLAQ EVLTPAQIKS ICQAILDSGK QYAIKKRKPFLMYSYYGTE YLGAAHGLSS ILQMLLSYHE HLKPSDRELV WQSVDFLMEQ EQNCNWPEL GETIERENEL VHWCHGAPGI AYLFAKAYLV SKKPQYLDTC IRCGELTWQK GLLKKGPGIC HGVAGSAYVF LLLYRLTGNS KYIYRAQRFA QFLFTEEFKA GSRVLESIYS LYEGFSGTVC FLIDLLQPNQ AEFPLFSVFV
Source	Mammalian cell
Target Names	LANCL3
Protein Names	Recommended name: LanC-like protein 3
Expression Region	1-420
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.