



Recombinant Human Pleckstrin homology domain-containing family M member 3 (PLEKHM3)

Product Code	CSB-MP741077HU
Abbreviation	PLEKHM3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ZWE6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEALEVDDIS PALEVTEEFF STLDNLEKA VQQAQEVYGIQ EVPELVGHEV LSNITDNGAM RNVTSLGKGG MIWDHCKSRL LETKAQNVFP AKEQFMVQRG TTPDNLSWME QKEASTFNFF NICQRRRDRP RSVNDLLDET STFKPGHARS RSDITQVDWR VVLKTTPLQQ QQQQQPLLQG PHVTRPSFLL PSPNKIEDAQ GNTEHKQTFP NILKKGYLEI RKDHDYSWQS CYAELSPYNL YFYSLDSSGN QNLATYQLS HFQISISVLGN LEARMVDTV L YDNTQLQLKA ESPWEALDWG QKLWEVVHAA VPGYMGRQNE LTISPGLGHH DDYTQNHFSQ KKTSGLLPPS PVL DSSKQYQ NILKSGTLYR LTVQNNWKAF TFVLSRAYLM AFQPGKLDED PLLSYNVDVC LAVQMDNLDG CDSCFQVIFP QDVLRLRAET RQRAQEWMEA LKIAANVARS SEQNLQVTLR NPKKDQMGGH ELRKNKRQSV TTSFLSILTT LSLERGLTAQ SFKCAGCQRS IGLSNGKAKV CNYSGWYYCS SCHVDDSFLI PARIVHNWDT SKYKVSKQAK EFLEYVYEEP LIDIQQENAM LYHHAEP LAA VLRLRQRLKS LRAYLFSCRA AVAEDLRRRI FPREYLLQQI HLYSLADLQQ VIEGKLAPFL GKVIKFATSH VYSCSLCSQK GFICEICNNG EILYPFEDIS TSRCESCGAV FHSECKEKS SV PCPRCVRREL QKKQKSFQWR LNMDESLEEA CTMFELSYQN T
Source	Mammalian cell
Target Names	PLEKHM3
Protein Names	Recommended name: Pleckstrin homology domain-containing family M member 3 Short name= PH domain-containing family M member 3
Expression Region	1-761
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.