



Recombinant Human Activator of 90 kDa heat shock protein ATPase homolog 2 (AHSA2)

Product Code	CSB-MP741128HU
Abbreviation	AHSA2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q719I0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAKWGQQGNPH WIVEEREDGT NVNNWRWTER DATSLSKGGK QELLVGIVVE NDAGRGEINE LKQVEGEASC SSRKGLIFF YEWNIKLGWK GIVKESGVKH KGLIEIPNLS EENEVDDTEV SLSKKKGDGV ILKDLMKTAG TAKVREALGD YLKALKTEFT TGMILPTKAM ATQELTVKRK LSGNTLQVQA SSPVALGVRI PTVALHMMEL FDTTVEQLYS IFTVKELTNK KIIMKWRCGN WPEEHYAMVA LNFVPTLGQT ELQLKEFLSI CKEENMKFCW QKQHFEEIKG SLQLTPLNG
Source	Mammalian cell
Target Names	AHSA2P
Protein Names	Recommended name: Activator of 90 kDa heat shock protein ATPase homolog 2
Expression Region	1-299
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.