



Recombinant *Desulfovibrio vulgaris* Exodeoxyribonuclease 7 small subunit (xseB)

Product Code	CSB-EP741239DDH-B
Abbreviation	xseB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q72CD5
Product Type	Recombinant Protein
Immunogen Species	<i>Desulfovibrio vulgaris</i> (strain Hildenborough / ATCC 29579 / DSM 644 / NCIMB 8303)
Purity	≥85% (SDS-PAGE)
Sequence	MAGRKKAADF EQQLARLQEI VDALEGGDLP LEKSVALYKE GLGLARASRE QLAKARNEIR LFTEGEVRDF DPEEGDDGDD R
Source	<i>E.coli</i>
Target Names	xseB
Protein Names	Recommended name: Exodeoxyribonuclease 7 small subunit EC= 3.1.11.6 Alternative name(s): Exodeoxyribonuclease VII small subunit Short name= Exonuclease VII small subunit
Expression Region	1-81
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.