



Recombinant *Lactobacillus johnsonii* Enolase 2 (eno2)

Product Code	CSB-BP741541LAW
Abbreviation	eno2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q74J64
Product Type	Recombinant Protein
Immunogen Species	<i>Lactobacillus johnsonii</i> (strain CNCM I-12250 / La1 / NCC 533)
Purity	≥85% (SDS-PAGE)
Sequence	MTVYVEKVRA LEIFDSRGNP TVEVHAYLSD GTVAKAEVPS GASTGEKEAV ELRDGGNRLQ GKGVTQAVTN VNGPINDALK GLSPYNQAEI DRTMIKLDGT LNKAKLGANA ILGTSMAIAR AAARSKDEPL YRYLGGCELE MPQTFHNVIN GGKHADNGID IQEFMITPVA KNSFRDGFKE IVNTYHALKA VIEEAGFETG LGDEGGFAPN LNSSEEALKM LRKAIKAGY KPRKDIAIAF DAAASSFYNT EDGKYHFEGH IWNGEEMLQY YDKLLKEFPE IISCEDPFDE NDWENFEKFT AKFGSTHQVV ADDNVCTNPK LVRKAIDKL CNSILIKLNQ IGTITETLET IRLARKNNMT TMVSHRSGET GDTFIADFTV ATNAGQLKSG APARSERVEK YNRLLEIENQ IGVENERLNH FPNNVDFD
Source	Baculovirus
Target Names	eno2
Protein Names	Recommended name: Enolase 2 EC= 4.2.1.11 Alternative name(s): 2-phospho-D-glycerate hydro-lyase 2 2-phosphoglycerate dehydratase 2
Expression Region	1-428
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.