



# Recombinant *Ashbya gossypii* AP-3 complex subunit mu (APM3)

<b>Product Code</b>	CSB-YP741585DOT
<b>Abbreviation</b>	APM3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q750L8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) ( <i>Eremothecium gossypii</i> )
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MYTSIYITDS KKNLVFEYLL TSQAPTFQQL CSKIAGRDAADAMVAISKD MSLYRQAVGP EKLWYWALCQ ACGDPLEPQM FLAQYHQVLI EYFDKEALTV KCLVNNADRL ALLHAMLDA GEVAVTDSNR LRQLVPLRND LSTILNSATK TLANTVKYAD SKQLFGAPVA TGKVEAGQTV PWRTADCRYV NNEIYVDLVE TVNATLRQKG SSLTLINGSL SGKIDVKCYL SGNPTVQLKL RTSGHPLDNS ALHRCVELGE AGVATMNFVP PDGRFTLAEY AIDLSAISQA ARRLTNLGLV TVSLASGLGQ HRDEFEIKVI IGNSTQVAEI EDLRITVYFP DISDAAKIKI LRTTHGGWES DLSRVRGVWA FDKQTAVGSV PVLRGCVENP ESTPHAPPVF PSHLAVSYSH VGQLPSGIRV DTIALSDLPP GSKPFKGVKY TSRAGDYIVR A
<b>Source</b>	Yeast
<b>Target Names</b>	APM3
<b>Protein Names</b>	Recommended name: AP-3 complex subunit mu Alternative name(s): AP-3 adapter complex mu3A subunit Adapter-related protein complex 3 mu subunit Mu-adaptin 3A Mu3-adaptin
<b>Expression Region</b>	1-451
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.