



Recombinant *Ashbya gossypii* RuvB-like helicase 1 (RVB1)

Product Code	CSB-EP741586DOT
Abbreviation	RVB1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q750R1
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MVQISEVKDQ VPGSSAGART AAHTHIKGLG LDEFGAAKQV EGGFVGQVEA REACGVIVDL IKAKRMSGRA ILLAGGPSTG KTALALAITQ ELGPKVPFCP LVGSELSVE VKKTETLMEN FRAIGLRIK EVKEVYEGEV TELTPEEAEN PLGGYGKTIS HVIVGLKSAK GTKTLRLDPT IYESIQREKV SVGDVIIYES NSGAVKRVGR SDAYATEFDL EAEYVPLPK GEVHKKKEII QDVTLHDLDV ANARPQGGQD VISMMGQLMK PKKTEITEKL RHEVNKVVAK YIDQGV AELV PGVLFIDEVN MLDIEIFTFL NRALELEIAP VVVLASNRGM TTVRGTEDVV SAHGIPDLI DRLLIVRTL P YTQDEIRVII EKRSKVENLQ LEQAALDLLA AMGSDMSLRY ALQLLTPAGI LAATAGRTEI LLSDIEEAKM LFLDAKRSTK ILESNSNYL
Source	E.coli
Target Names	RVB1
Protein Names	Recommended name: RuvB-like helicase 1 EC= 3.6.4.12
Expression Region	1-459
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.