



# Recombinant *Ashbya gossypii* Protein ESC2 (ESC2)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-YP741620DOT   |
| <b>Abbreviation</b>      | ESC2  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q756K3  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | <i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) ( <i>Eremothecium gossypii</i> )  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MGEFNEYDDI DDFFCNDISP LHDEADEFSG TDIALPSQLY GKGPKPRAGE<br>SKALRDRNKE SGSGETEAK GDTLATKRGR SSWRRGASES SRSSSLGSSS<br>PSDSSSGRSL SPVRPAKRKR TDEQQPESEA NRFLAEMERD IELSAGPGGG<br>ETATGKGGRD TDARVYNVGF ISRIEGSRNR RVNVKVTGQK QFATFLAIAL<br>AAFSKQHNIR KSLKPKYQAD QVKIYREGVE VFKFMTCDSE NIEEPYNASA<br>TDIEVYIVPV DEATAFEQEW KRKFEERVLM LSNSSFLEVM DDIEDDNDDF<br>LVNEYEKALV NARSLNETEL AVREGTPADE SSQLLKIVLL GSDNKKVFIH<br>VRPTTLLKV AEHYRVAKEL PPTVQLSLMF DHEEIDLDDT ICNIDIEDGD<br>IIEVVVK |
| <b>Source</b>            | Yeast   |
| <b>Target Names</b>      | ESC2  |
| <b>Protein Names</b>     | Recommended name: Protein ESC2 Alternative name(s): Establishes silent chromatin protein 2  |
| <b>Expression Region</b> | 1-407   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life   |



of lyophilized form is 12 months at -20°C/-80°C.