



Recombinant *Ashbya gossypii* Phosphatidylinositol transfer protein CSR1 (CSR1)

Product Code	CSB-YP741625DOT
Abbreviation	CSR1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q757H2
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MTQTAAPGRV QTLTKEQEKV LKQVWVHLFQ FWGIDVDGSA VLTAKEPEPA APSGKGRKLL GLFGKRKDAT GAARKESNGK AAKVYDAEK VEDSDAEKEK PTPQKVEGLE EMYELLKELD GAAVSKEFWS MLRCDYPDNL LLRFVRARKW DINKAMIMMA HSLRWRLNEG KPEDIVFGGE RGAQKADKKG IVKQLELGKA TVRGFDKNGC PIVYVRPRLH HAADQTEAET SEYSLIIIEQ ARLFLKEPCD TATILFDLSG FSMANMDYAP VKFLITCFEA HYPECLGKLF IHKAPWIFPP IWNIIKNWLD PVVAAKIAFT KTAADLEEFI PAEQIPLELG GKDEYNFDGF VMPDGSADTK LSDEKGKAAV LEEREAIIKR FIDATISWIE STSDEESAKW LEKKIDLSKE LSENYSKLDP YIRSRSFYDV NGTLKV
Source	Yeast
Target Names	CSR1
Protein Names	Recommended name: Phosphatidylinositol transfer protein CSR1
Expression Region	1-436
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.