



Recombinant *Ashbya gossypii* 54S ribosomal protein L2, mitochondrial (MRPL2)

Product Code	CSB-BP741680DOT
Abbreviation	MRPL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q75ET6
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	ATKRAA GSRTSMKDSA GRRLGPKKTE GQRVEVGQII MRQRGTRFYP GENVGIGKDH TLFALEPGWV RYYLDPFHTG RKFVGVVALYQ DMRLPIEHFA PRVRRFGRQL LSGKEAETEE QALPRSVFLA KEKILERAQQ RADAREQRRR EFGHILREEL GLALDQDAEQ LATEYLLRVQ TNLKSGFNSA DARFNAMYHM EVRMRNTPEK VGQMDLLKKT VASVDAATSF SNKFELGRHI SEGERVAWRE ALHKDLAGLV IRTAEEKQRV VERLKEASKY LSLSEEIHLR RKFLKPVKPE TEAVAEAAGK GTVTIKRFNY GTGKVDTIIR EKKAFLAKL
Source	Baculovirus
Target Names	MRPL2
Protein Names	Recommended name: 54S ribosomal protein L2, mitochondrial
Expression Region	25-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.