



Recombinant *Oryza sativa* subsp. japonica Protein DROOPING LEAF (DL)

Product Code	CSB-MP741726OFG
Abbreviation	DL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q76EJ0
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MDLVSPSEHL CYVRCTYCNT VLA VGV PCKR LMDTVTVKCG HCN NLSFLSP RPPMVQPLSP TDHPLGPFQG PCTDCRRNQP LPLVSPTSNE GSPRAPFVVK PPEKKHRLPS AYNRFMREEI QRIKAAKPGI PHREAFSMAA KNWAKCDPRC SSTVSTNSN PEPRVVA API PHQERANEQV VESFDIFKQM ERSG
Source	Mammalian cell
Target Names	DL
Protein Names	Recommended name: Protein DROOPING LEAF Alternative name(s): Protein CRABS CLAW homolog Short name= Protein CRC homolog
Expression Region	1-194
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.