



Recombinant Newcastle disease virus Nucleoprotein (N)

Product Code	CSB-EP741755NCT-B
Abbreviation	N
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q77K03
Product Type	Recombinant Protein
Immunogen Species	Newcastle disease virus (strain Chicken/United States/B1/48) (NDV)
Purity	>85% (SDS-PAGE)
Sequence	MSSVFDEYEQ LLAAQTRPNG AHGGGEKGST LKVDVPVFTL NSDDPEDRWS FVVFCLRIAV SEDANKPLRQ GALISLLCSH SQVMRNHVAL AGKQNEATLA VLEIDGFANG TPQFNNSGV SEERAQRFAM IAGSLPRACS NGTPFVTAGA EDDAPEDITD TLERILSIQA QVWVTVAKAM TAYETADESE TRRINKYMQQ GRVQKKYILY PVCIRSTIQLT IRQSLAVRIF LVSELKRGRN TAGGTSTYYN LVGDVDSYIR NTGLTAFFLT LKYGINTKTS ALALSSLSGD IQKMKQLMRL YRMKGDNAPY MTLGSDSDQM SFAPAEYAQL YSFAMGMASV LDKGTGKYQF AKDFMSTSFV RLGVEYAQAQ GSSINEDMAA ELKLTPAARR GLAAAAQRVS EVTSSIDMPT QQVGVLTGLS EGGSQALQGG SNRSQQQPEA GDGETQFLDL MRAVANSMRE APNSAQGTPQ SGPPPTPGPS QDNDTDWGY
Source	E.coli
Target Names	N
Protein Names	Recommended name: Nucleoprotein Alternative name(s): Nucleocapsid protein Short name= NP Short name= Protein N
Expression Region	1-489
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.