



# Recombinant Newcastle disease virus Nucleoprotein (N)

<b>Product Code</b>	CSB-YP741755NCT
<b>Abbreviation</b>	N
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q77K03
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Newcastle disease virus (strain Chicken/United States/B1/48) (NDV)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSSVFDEYEQ LLAAQTRPNG AHGGGEKGST LKVDVPVFTL NSDDPEDRWS FVVFCLRIAV SEDANKPLRQ GALISLLCSH SQVMRNHVAL AGKQNEATLA VLEIDGFANG TPQFNNSGV SEERAQRFAM IAGSLPRACS NGTPFVTAGA EDDAPEDITD TLERILSIQA QVWVTVAKAM TAYETADESE TRRINKYMQQ GRVQKKYILY PVCIRSTIQLT IRQSLAVRIF LVSELKRGRN TAGGTSTYYN LVGDVDSYIR NTGLTAFFLT LKYGINTKTS ALALSSLSGD IQKMKQLMRL YRMKGDNAPY MTLGSDSDQM SFAPAEYAQL YSFAMGMASV LDKGTGKYQF AKDFMSTSFV RLGVEYAQAQ GSSINEDMAA ELKLTPAARR GLAAAQRVS EVTSSIDMPT QQVGVLTGLS EGGSQALQGG SNRSQQQPEA GDGETQFLDL MRAVANSMRE APNSAQGTPQ SGPPPTPGPS QDNDTDWGY
<b>Source</b>	Yeast
<b>Target Names</b>	N
<b>Protein Names</b>	Recommended name: Nucleoprotein Alternative name(s): Nucleocapsid protein Short name= NP Short name= Protein N
<b>Expression Region</b>	1-489
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.