



Recombinant Escherichia coli Putative uncharacterized protein yghX (yghX)

Product Code	CSB-EP741907ENV
Abbreviation	yghX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7DFU6
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	MTITLTGKSR NVSSSTLRRS AVGGMTALAL FDLLKPNYAL ATQVEFTDPE IVA EYITYPS PNGHGEVRGY LVKPAKMSGK TPAVVVVHEN RGLNPYIEDV ARRVAKAGYI ALAPDGLSSV GGYPGNDDKG RELQQQVDPT NLMNDFFAAI EFMQRYPQAT GKVGITGFCY GGGVSNAAAV AYPELACAVP FYGRQAPTAD VAKIEAPLLL HFAELDTRIN EGWPAYEAAL KANNKVYEAY IYPGVNHGFH NDSTPRYDKS AADLAWQRTL KWFDKYLS
Source	E.coli
Target Names	yghX
Protein Names	Recommended name: Putative uncharacterized protein yghX
Expression Region	1-278
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.