



Recombinant *Oryza sativa* subsp. japonica Probable E3 ubiquitin-protein ligase XBOS35 (XBOS35)

Product Code	CSB-EP741910OFG-B
Abbreviation	XBOS35
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7EZ44
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MGLLGMVGDG FGCSATGERL VSAARDGDIQ EAMALLELNP RLARYSTFGI RNSPLHYSAA KGHHEIVSLL IESGVDINLR NCRGQTALMQ ACLYGHWKVV QILVLFKANI HKKDCFSGAT AIHFAALKGH TRCLRLLVAD YVPSLPEFWS VMHAKCTDET NKEAFDAVAL RRLINNKSDG GVTPLHLAAL HGHAECVQLL LDLGASVSEV TINDGSTIDL IGSGSTPLHY AACGGSVAVCC QLLVAAGANM RAQNTNGLTP LMVARSWHKS SVEGILTKRS EVPVRILPSS YLSLPLMSIV KIARECGWRK TSVSSVCHDP CAICLDTECT VSAEGCGHEF CTKCALYLCA TASSSTSIRG VPGSIPCLC RHTIVSFVRL TSTTPIKELP WTNKSLALCA AGASTGSKYA GPAAITSSKY AGSLHRRSEM RSLRSSVVDL GCSSFRTASS GKLSSIKLNC TGADETMPCL VNCFRPDVQR SSSYRERIRR YSQFS
Source	E.coli
Target Names	XBOS35
Protein Names	Recommended name: Probable E3 ubiquitin-protein ligase XBOS35 EC= 6.3.2.- Alternative name(s): Ankyrin repeat domain and RING finger-containing protein XBOS35 XB3 protein homolog 5
Expression Region	1-495
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.