



# Recombinant *Oryza sativa* subsp. japonica Probable E3 ubiquitin-protein ligase XBOS35 (XBOS35)

<b>Product Code</b>	CSB-MP7419100FG
<b>Abbreviation</b>	XBOS35
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7EZ44
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGLLGMVGDG FGCSATGERL VSAARDGDIQ EAMALLELNP RLARYSTFGI RNSPLHYSAA KGHHEIVSLL IESGVDINLR NCRGQTALMQ ACLYGHWKVV QILVLFKANI HKKDCFSGAT AIHFAALKGH TRCLRLLVAD YVPSLPEFWS VMHAKCTDET NKEAFDAVAL RRLINNKSDG GVTPLHLAAL HGHAECVQLL LDLGASVSEV TINDGSTIDL IGSGSTPLHY AACGGSVAVCC QLLVAAGANM RAQNTNGLTP LMVARSWHKS SVEGILTKRS EVPVRILPSS YLSLPLMSIV KIARECGWRK TSVSSVCHDP CAICLDTECT VSAEGCGHEF CTKCALYLCA TASSSTSIRG VPGSIPCLC RHTIVSFVRL TSTTPIKELP WTNKSLALCA AGASTGSKYA GPAAITSSKY AGSLHRRSEM RSLRSSVVDL GCSSFRTASS GKLSSIKLNC TGADETMPCL VNCFRPDVQR SSSYRERIRR YSQFS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	XBOS35
<b>Protein Names</b>	Recommended name: Probable E3 ubiquitin-protein ligase XBOS35 EC= 6.3.2.- Alternative name(s): Ankyrin repeat domain and RING finger-containing protein XBOS35 XB3 protein homolog 5
<b>Expression Region</b>	1-495
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.