



Recombinant *Oryza sativa* subsp. japonica Peptide-N (4)- (N-acetyl-beta- glucosaminyl)asparagine amidase (PNG1)

Product Code	CSB-BP741911OFG
Abbreviation	PNG1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7F0R1
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MVARRFVVRQ GGGGGGGGEA EEHEVEYDTE HGLDILRLQI FSLTSVPPEL QKIVVEADGS VVDDGTDLEA ISEGLRLVAI TGEAAAAEAA AAAEAARAQE KSDEELARMI QAEEEEALLQ QYSIRNDGGE EFRERVEPYM HQVLMYEDPM RQEAARKTVP MDELQEKALV SLAKEGNFSP SKDEEDHAFL LQLLFWFKQS FRWVNAPPCD SCGRETFNVG MGTALPSEIK FGANRVEIYR CNYCSSTTRF PRYNDPYKLL ETRKGRCGEW ANCFTFYCRS FGYEARLILD FTDHVVTECF SNLYGRWMHL DPCEGVYDNP LLYEKGWNKK LDYVIAISKD GVRDVTKRYT RKWHEVLSRR IITSEDTVSA ILSSITGKYR SGLSIDGLTA LENRDKKESE ELSKAAYLEV DTSISLPGRQ SGSVEWRKAS QKCSTYILSI TSGNGCG
Source	Baculovirus
Target Names	PNG1
Protein Names	Recommended name: Peptide-N(4)-(N-acetyl-beta-glucosaminyl)asparagine amidase EC= 3.5.1.52 Alternative name(s): Peptide:N-glycanase
Expression Region	1-447
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.