



# Recombinant *Oryza sativa* subsp. japonica Peptide-N (4)- (N-acetyl-beta- glucosaminyl)asparagine amidase (PNG1)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-YP741911OFG   |
| <b>Abbreviation</b>      | PNG1  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q7F0R1  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | <i>Oryza sativa</i> subsp. japonica (Rice)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MVARRFVVRQ GGGGGGGGEA EEHEVEYDTE HGLDILRLQI FSLTSVPPEL<br>QKIVVEADGS VVDDGTDLEA ISEGLRLVAI TGEAAAAEAA AAAEAARAQE<br>KSDEELARMI QAEEEEALLQ QYSIRNDGGE EFRERVEPYM HQVLMYEDPM<br>RQEAARKTVP MDELQEKALV SLAKEGNFSP SKDEEDHAFL LQLLFWFKQS<br>FRWVNAPPCD SCGRETFNVG MGTALPSEIK FGANRVEIYR CNYCSSTTRF<br>PRYNDPYKLL ETRKGRCGEW ANCFTFYCRS FGYEARLILD FTDHVVTECF<br>SNLYGRWMHL DPCEGVYDNP LLYEKGWNKK LDYVIAISKD GVRDVTKRYT<br>RKWHEVLSRR IITSEDTVSA ILSSITGKYR SGLSIDGLTA LENRDKKESE<br>ELSKAAYLEV DTSISLPGRQ SGSVEWRKAS QKCSTYILSI TSGNGCG |
| <b>Source</b>            | Yeast   |
| <b>Target Names</b>      | PNG1  |
| <b>Protein Names</b>     | Recommended name: Peptide-N(4)-(N-acetyl-beta-glucosaminyl)asparagine amidase EC= 3.5.1.52 Alternative name(s): Peptide:N-glycanase   |
| <b>Expression Region</b> | 1-447   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.  |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.