



Recombinant *Oryza sativa* subsp. *japonica* Probable NAD (P)H-dependent oxidoreductase 2 (Os10g0113100, LOC_Os10g02380, LOC_Os10g02390)

Product Code	CSB-EP741918OFG
Abbreviation	Os10g0113100, LOC_Os10g02380, LOC_Os10g02390
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7G765
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MAMATIPEVP ASALLPTMPR IGMGTAAFPF TSSEETTAAL LRAIELGYRH FDTARLYATE GCVSEAVAEA VRRGLVASRA DVFVTSKLWC SDLHAGRVP AARETLRNLG MDYVDLLLH WPATVAPGSY DFPFPKEEMA PAFDMEGVWR GMEECHRLGL ARAIGVSNFS AKKLEQLLSF AVVRPAANQV EMNPMWQQRT LREVCRRREGV QLCGYSPLGA KGTPWGSAAV MDSGVLHDIA QTKGKTLAQI CLRWMYEQGD VLLVKTYNEN RMKENLDIFD WELTEEERDK ISKLPQQRGL TGMQFVCDNG PYKCVEDLWD GA
Source	E.coli
Target Names	Os10g0113100
Protein Names	Recommended name: Probable NAD(P)H-dependent oxidoreductase 2 EC= 1.- .-.-
Expression Region	1-322
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.