



Recombinant Drosophila melanogaster Vacuolar protein-sorting-associated protein 25 (Vps25)

Product Code	CSB-EP741957DLU-B
Abbreviation	Vps25
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7JXV9
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MAEFQWPWEY TFPPFFTLQP HEETRQQQLK VWTDLFLKYL RHTNRFTLSI GDQNSPLFHN EALKRRLSPE LVLAILGELE RSGHANPLDK RRQEWQVYWF TLEEYGNMVF DWVQETGQTN TICTLYEIAS GENTSHLDFY GVDEAVLLSA LRLLEEKGRG ELIEMDGSHG VKFF
Source	E.coli
Target Names	Vps25
Protein Names	Recommended name: Vacuolar protein-sorting-associated protein 25 Alternative name(s): ESCRT-II complex subunit VPS25 Vacuolar protein sorting 25
Expression Region	1-174
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.