



Recombinant Human BTB/POZ domain-containing protein KCTD9 (KCTD9)

Product Code	CSB-BP741970HU
Abbreviation	KCTD9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7L273
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MRRVTFLFLNG SPKNGKVVAV YGTLSDLLSV ASSKLGIKAT SVYNGKGGLI DDIALIRDDD VLFVCEGEPF IDPQTDSKPP ELLGFHTDW LTLNVGGRYF TTTRSTLVNK EPDSMLAHMF KDKGVWGNKQ DHRGAFLIDR SPEYFEPILN YLRHGQLIVN DGINLLGVLE EARFFGIDSL IEHLEVAIKN SQPPEDHSPI SRKEFVRFL ATPTKSELRC QGLNFSGADL SRLDLRYINF KMANLSRCNL AHANLCCANL ERADLSGSVL DCANLQGVKM LCSNAEGASL KLCNFEDPSG LKANLEGANL KGVDMEGSQM TGINLRVATL KNAKLKNCNL RGATLAGTDL ENCDLSGCDL QEANLRGSNV KGAIFFEMLT PLHMSQSVR
Source	Baculovirus
Target Names	KCTD9
Protein Names	Recommended name: BTB/POZ domain-containing protein KCTD9
Expression Region	1-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.