



Recombinant Human G2/M phase-specific E3 ubiquitin-protein ligase (G2E3)

Product Code	CSB-EP741972HU-B
Abbreviation	G2E3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7L622
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MNESKPGDSQ NLACVFCRKH DDCPNKYGEK KTKEKWNLTV HYCLLMSSG IWQRGKEEEG VYGFLIEDIR KEVNRASKLK CCVCKKNGAS IGCVAPRCKR SYHFPCGLQR ECIFQFTGNF ASFCWDHRPV QIITSNNYRE SLPCTICLEF IEPISYNIL RSPCKNAWF HRDCLQVQAI NAGVFFFRCT ICNNSDIFQK EMLRMGIHIP EKDASWELEE NAYQELLQHY ERCDVRRRCRC KEGRDYNAPD SKWEIKRCQC CGSSGTHLAC SSLRSWEQNW ECLECRGIY NSGEFQKAKK HVL PNSNVG ITDCLLEESS PKLPRQSPGS QSKDLLRQGS KFRRNVSTLL IELGFQIKKK TKRLYINKAN IWNSALDAFR NRNFNPSYAI EVAYVIENDN FGSEHPGSKQ EFLSLLMQHL ENSSLFEGSL SKNLSLNSQA LKENLYYEAG KMLAISLVHG GPSPGFFSKT LFNCLVYGPE NTQPILDDVS DFDVAQIIIR INTATTVADL KSIINECYN YLELIGCLRLI TTLSDKYMLV KDILGYHVIQ RVHTPFESFK QGLKTLGVLE KIQAYPEAFC SILCHKPESL SAKILSELFVHTLPDVKAL GFWNSYLQAV EDGKSTTTME DILIFATGCS SIPPAGFKPT PSIECLHVDF PVGNKCNNCL AIPITNTYKE FQENMDFTIR NLRLEKEES SHYIGH
Source	E.coli
Target Names	G2E3
Protein Names	Recommended name: G2/M phase-specific E3 ubiquitin-protein ligase EC=6.3.2.-
Expression Region	1-706
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.