



Recombinant Human CXXC-type zinc finger protein 5 (CXXC5)

Product Code	CSB-BP741975HU
Abbreviation	CXXC5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7LFL8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSSLGGGSQD AGGSSSSSTN GSGGSGSSGP KAGAADKSAV VAAAAPASVA DDTPPPERRN KSGIISEPLN KSLRRSRPLS HYSSFGSSGG SGGGSMMGGE SADKATAAAA AASLLANGHD LAAAMAVIDKS NPTSKHKSGA VASLLSKAER ATELAAEGQL TLQQFAQSTE MLKRVVQEHL PLMSEAGAGL PDMEAVAGAE ALNGQSDFPY LGAFPINPGL FIMTPAGVFL AESALHMAGL AEYPMQGELA SAISSGKKKR KRCGMCAPCR RRINCEQCSS CRNRKTGHQI CKFRKCEELK KKPSAALEKV MLPTGAAFRW FQ
Source	Baculovirus
Target Names	CXXC5
Protein Names	Recommended name: CXXC-type zinc finger protein 5 Alternative name(s): Putative MAPK-activating protein PM08 Putative NF-kappa-B-activating protein 102
Expression Region	1-322
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.