



Recombinant Chalinolobus morio Hemoglobin subunit alpha (HBA)

Product Code	CSB-MP741989CHM
Abbreviation	HBA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7M2Y4
Product Type	Recombinant Protein
Immunogen Species	Chalinolobus morio (Chocolate-wattled bat) (Lobe-lipped bat)
Purity	>85% (SDS-PAGE)
Sequence	VLSPADKSNV KAAWDKVGGS AGDYGAEALE RMFLSFPTTK TYFPHFDLSH GSAQVKGHGK KVGDALGNAV AHLDDLPGAL SALSDLHAHK LRVDPVNFKL LSHCLLVTV ACHHPNDFTPA VHASLDKFLA NVSTVLVSKY R
Source	Mammalian cell
Target Names	HBA
Protein Names	Recommended name: Hemoglobin subunit alpha Alternative name(s): Alpha-globin Hemoglobin alpha chain
Expression Region	1-141
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.