



Recombinant Notophthalmus viridescens Fibroblast growth factor 1 (fgf1)

Product Code	CSB-EP742452NEO-B
Abbreviation	fgf1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7SIF8
Product Type	Recombinant Protein
Immunogen Species	Notophthalmus viridescens (Eastern newt) (Triturus viridescens)
Purity	>85% (SDS-PAGE)
Sequence	QKPKLLYCSN GGYFLRIFPD GKVDGTRDRS DPYIQLQFYA ESVGEVYIKS LETGQYLAMD SDGQLYASQS PSEECFLER LEENNYNTYK SKVHADKDFW VGIKKNKGTK PGSRTHFGQK AILFLPLPVS SD
Source	E.coli
Target Names	fgf1
Protein Names	Recommended name: Fibroblast growth factor 1 Short name= FGF-1 Alternative name(s): Acidic fibroblast growth factor Short name= aFGF Heparin-binding growth factor 1 Short name= HBGF-1
Expression Region	1-132
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.